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## Technical specifications

	Rango presión de entrada Inlet pressure range	Máxima presión de salida Maximum pressure outlet	Caudal Flow	Peso Weight
OXYMO P-15 MICRO	2.8 – 7 bar	1.5 ± 0.15 bar	> 25 L/min	150 g
OXYMO P-15	2.8 – 7 bar	1.5 ± 0.15 bar	> 25 L/min	300 g


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

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## Oxygen regulators OXYMO P15 and P15Micro Directions for use

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OXYMO P-15 Micro and OXYMO P-15 are accessories to the device OXYMO. This manual should be read in conjunction with the OXYMO Manual Ref. 9550600

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### WARRANTY

HERSILL warrants to the purchaser that this device is free from defects in material and workmanship for a period of 2 YEARS from the date of purchase by the original user. The guarantee does not cover damage to the device due to improper use or use of parts unauthorised by HERSILL. This warranty does not cover normal wear and tear, staining, discoloration or other cosmetic irregularity which does not impede or degrade the normal product performance

HERSILL, S.L.  
Puerto de Navacerrada 3 - P.I. Las Nieves  
28935 Móstoles (Madrid) – Spain  
[www.hersill.com](http://www.hersill.com)

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## 1. Purpose. Intended use

OXYMO P-15 and P-15 MICRO, are low pressure regulators for medical oxygen, and should be exclusively used with the demand oxygen delivery system OXYMO

**CONTRAINDICATIONS** : None known

## 2. Cautions and warnings

-  Use this equipment only for the purpose for which it was designed
-  Do not smoke, and stay away from any flame or heat source, while using the device or you are close to someone who uses it.
-  Keep away from the device every substance that can be flammable, as lubricants, petroleum products, alcohols, aerosols, etc
-  Prevent the device can be manipulated by unauthorized or untrained persons
-  Federal laws in USA restricts this device to sale by or on the order of a physician (CFR 21 801.109.
-  Disposal of device, components and accessories must strictly conform to local laws and regulations.

## 3. Operating instructions (Fig.1 y Fig. 2)

For the oxygen supply what OXYMO need is right, follow the instructions below:

- Connect the equipment to a source of pressure between 2.8 and 7 bar, with the inlet (2) or (11), as appropriate
  -  The valve of the pressure source and the inlet of the regulator must be of the same type. Consult the different types available (see Fig. 5)
  -  The regulator must be connected to a pressure source. Be sure **NOT** to connect to an outflow
- Introducing connector (5) of the spiral hose (4) provided with the OXYMO equipment (7), the outlet valve (3) or (12) of the regulator, as appropriate.
- Connecting the valve (6) of the hose, to the inlet (8) of OXYMO (7)

### \* Note 1.

In the case of the OXYMO P-15 regulator (1), you can hang it on OXYMO (7) through the support (9) (see Fig. 3)

### \* Note 2.

In case of failure of the OXYMO (7), you can use the regulators OXYMO P-15 (1) and OXYMO P-15 MICRO (10) for supplying a continuous flow of 3 L/min with the adapter (13) (See Fig. 4), which is connected directly to the nasal cannula.

### \* Note 3

In the OXYMO, if the alarm for low supply pressure is activated (see chap. 8.2, OXYMO user manual), verify that there is not a fault in the oxygen supply, and that the connection between regulator, the pressure source and the OXYMO is correct.

If the problem persists, perform the function test (see below)

## 4. Function test

- a) Check with a manometer, if the pressure in the regulator outlet, when the flow is "0", it agree with that indicated in the table of technical specifications.
- b) Check that the maximum flow is able to give the regulator is higher than indicated in the table of technical specifications.

## 5. Maintenance

- The equipment does not require a specific maintenance
- It is recommended to check the device at least once every five years.
- Repairs must be performed by HERSILL or authorized personnel.